

## . Introduction

BREATH

BIOPSY

Exhaled breath contains hundreds of volatile organic compounds (VOCs), including those originating from biological processes in the body, and therefore have the potential to serve as biomarkers for clinical applications.

To reduce the variability and challenges associated with biomarker development, utilizing mice in controlled laboratory settings is necessary for expediting the identification and validation of breath biomarkers for clinical use.

Therefore, we have developed a reliable system for studying mouse breath in the laboratory. The platform we provide enhances on-breath VOC signals by significantly reducing background noise. This will improve foundational studies of breath and expedite the translation of breath analysis to the clinic.

### 2. Methods

A total of 15 intubated C57BL/6JRj mice breath samples were collected using a modified flexiVent<sup>®</sup> small animal ventilator (SCIREQ, Montreal, Quebec, Canada). An ambient filter and a flexiVent filter were both connected to the ventilator by flexible tubing and a sorbent tube for breath collection (Figure 1).

To compare mice and human breath profiles, a total of 13 healthy human volunteers were recruited. Each study participant provided a single breath sample collected using the ReCIVA® Breath Sampler.

For both mice and human samples, system blanks were collected to help distinguish breath from background compounds. All samples were analyzed using the Breath Biopsy® OMNI® standardized breath analysis method via thermal desorption gas chromatography and mass spectrometry (TD-GC-MS).

VOCs were identified through in-house HRAM and external NIST libraries, with compound IDs assigned following the Metabolomics Standard Initiative (MSI) standards.

Three metrics- standard deviation, paired t-tests and receiver operating characteristic area under the curve (ROC-AUC) were used in combination to determine VOCs 'on-breath' from background.



# **Development of a new breath collection method for analyzing volatile organic** compounds from intubated mouse models

### **3. Results**



*Figure 2 -* The different levels of total relative VOC concentration in different blank samples to test the background of the flexiVent system. The control was an empty tube, A was ambient air in the room, B was the samples. filtered ambient air, C was the air in the flexiVent system, and D was the filtered air in the flexiVent system.

To optimize untargeted analysis, a clean background is necessary to maximize the signal-to-noise ratio. Results (Figure 2) show that a second filter removes flexiVent contamination, identifying signals of interest more easily.

When comparing background signal between human blanks and mice blanks, 498 out of 661 (75.34%) features were significantly different. Most of the human blanks have higher signal than mice blanks, demonstrating flexiVent system have sufficiently lower background contamination levels for untargeted analysis of mouse breath (Figure 3).



*Figure 4 -* The number of overlapping identified "on-breath" VOCs by each metric type.

Depending on which classification thresholds were used, 16 – 73 on-breath VOCs were identified in mice (Figure 4). To assess translatability, we compared and found 57 common on-breath VOCs between mouse and human (Figure 5). While more VOCs were classified as on-breath in humans, this was expected due to the differences in relative size, lung volume, and environmental exposure between mice and humans.

Figure 1 - A

schematic showing the mouse breath sampling system used for this study. This system consists of a flexiVent small animal ventilator, an ambient filter, and a flexiVent filter both connected to the ventilator by flexible tubing, and a sorbent tube for breath collection.

Alastair Taylor<sup>1</sup>, Sylvia Blum<sup>2</sup>, Madeleine Ball<sup>1</sup>, Owen Birch<sup>1</sup>, Hsuan Chou<sup>1</sup>, Julia Greenwood<sup>1</sup>, Shane Swann<sup>1</sup>, Lara Pocock<sup>1</sup>, Max Allsworth<sup>1</sup>, Billy Boyle<sup>1</sup>, Kerstin Geillinger-Kaestle<sup>2</sup>



Figure 3 - Volcano plot comparing all VOCs in human and mouse blank

on-breath compounds identified in mouse breath and human breath by fold change compared to blank concentration. The off-breath compounds are shown in grey, those only in mice in orange, those only in humans in purple, and the on-breath compounds in both mice and humans are shown in blue.

Figure 5 - A plot of the



*Figure 6 - Examples of identified compounds. Trimethylamine (TMA) and dimethyl sulfone are identified* on-breath in both humans and mice. 2-Butanol and methyl nitrate are on-breath in mice only.

Many of the compounds common to mouse and human breath were associated with the gut microbiome, including TMA and dimethyl sulfone. Some compounds appear to be exclusively on-breath in mice only, including methyl nitrate and 2-butanol, both with a high signal intensity. These VOCs rarely observed in humans may represent a species-specific difference in breath composition or reflect a difference in environmental or dietary exposure.

### **4.** Conclusions

The results of this study present a reliable mouse breath sampling and analysis platform that can be used to compare the composition of mouse breath with human breath and establishes mice as a viable animal model for the pre-clinical study of breath biomarkers.

### owlstonemedical.com

<sup>1</sup>Owlstone Medical Ltd., Cambridge, Cambridgeshire, UK <sup>2</sup>Boehringer Ingelheim, Biberach, Germany \*email: breathbiopsy@owlstone.co.uk

• We identified 472 compounds in mouse breath, and 73 (15.47%) of these were considered as 'on-breath'. This demonstrates our three quantitative metrics are capable to distinguish signals that suggest "on-breath" compounds from background contamination.

• When comparing on-breath compounds, 57 (29.08%) were common between mouse and human, and some were linked to suspected biological functions.