

Breath-based biomarker detection

informing drug development and future treatment regimes

The recent approval of the first liquid biopsy test for the diagnosis of non-small cell lung cancer (NSCLC) demonstrates how biomarker-based detection tools are becoming important components of precision medicine-based drug treatment regimes.

Precision medicine allows detailed insights into the disease processes that underlie clinical disease. This field has first started with genetics, and genetic mutation-based models that can predict disease state such as tumour phenotype or likely therapeutic response are gaining much ground in oncology.

A genetic focus to precision medicine is limited, however, as this is a largely static picture that does not capture the biological variability of disease. In recent years there has been increasing focus on downstream molecules through transcriptomics, proteomics and metabolomics. Metabolomics is particularly interesting to stratify patients to the right treatment and inform development pipelines in the future, as analysis can be performed on molecules closely linked to disease activity.

Furthermore, the administration of drugs results in additional metabolites which can be used to evaluate pharmacokinetics and therapy efficacy. This paper reviews the use of the volatile fraction of

metabolites, so-called volatile organic compounds (VOCs), as a means to characterise disease processes and link the right patient to the right treatment.

Biomarkers for disease: example of liquid biopsy

Novel cancer assessment techniques such as liquid biopsy are becoming increasingly common. Liquid biopsy is a test done on a sample of blood to look for cancer cells from a tumour that are circulating in the blood, or for pieces of DNA from tumour cells that are in the blood. Results from this type of biomarker test can be returned to the patient in a matter of days, greatly reducing delay and uncertainty from traditional biopsies. Replacing invasive procedures with simple blood tests will also allow clinicians to greatly increase compliance rates. The applications include selection of therapies matched to individual patient mutations as well as screening and early detection.

While liquid biopsy is undoubtedly a useful tool

**By Billy Boyle and
Dr Mike Murphy**

there is still little agreement in identified breath biomarkers within a given disease, and breath testing is rarely used routinely in clinical settings.


New improvements in our ability to collect and inexpensively analyse breath opens up a range of new applications for breath biomarker tests. Improvements combining high sensitivity to VOCs and flexibility and reproducibility in sample collection allow breath biomarkers to be incorporated into early clinical trials and also to be used to screen and stratify patients for those new therapies.

The traditional breath collection technique of exhaling directly on to a sensor risks swamping the VOC signal with humidity and volatile chemicals naturally present in the background environment. Breath storage is most commonly carried out by exhaling into Tedlar® bags, which suffer from short breath sample retention times and poor transportability, and are not always suitable for patients with impaired breathing. To address these issues, a consortium of breath scientists (<http://www.breathe-free.org/>) developed a breath sampler that enables quick and easy breath collection in patients. This sampler allows robust, targeted selection of different portions of breath, allowing a focus on those fractions carrying the highest concentration of target biomarkers. Combining such a sampler with a clean air supply eliminates humidity and background contamination issue, and allows for comfortable sample collection.

VOCs and translational medicine

With improved sample collection, it now becomes feasible to introduce breath biopsy into the clinical setting. To perform the analysis on the captured breath sample, there are already high end analytical instruments such as gas chromatograph mass spectrometers (GC-MS) that are capable of detecting VOCs with high analytical accuracy. A GC-MS is the gold standard in sample identification, where chemicals in the sample are first separated by their chemical properties and then ionised and further distinguished by their mass to charge ratio. While their utility is well-proven, the cost and complexity of these machines restrict translation deployment in the clinic and make them an unlikely solution in the primary care setting.

On the other end of the scale, there are a plethora of simple materials-based sensors (so-called 'Electronic noses', resonant cantilevers, nanotubes, and other devices) that are cost-effective and easy to use. They operate by pattern recognition: this allows them to differentiate between 'disease' and 'control' states but does not give the identity of



Breath allows us to sample large volumes of blood completely non-invasively, leading to improved test sensitivity


the specific biomarkers that drive the detection signal, and so cannot help identify the mechanism behind the disease state.

A Field Asymmetric Ion Mobility Spectrometry (FAIMS) sensor platform combines the analytical power of a GC-MS with the flexibility and ease-of-use of materials-based sensors. The FAIMS chip can be programmed in software to detect targeted disease VOCs with high sensitivity and selectivity. First, the VOCs are ionised allowing their path to be altered by electric fields. The second stage of the process takes advantage of this by applying an alternating voltage across the channels of the FAIMS chip that the ions are travelling through. This creates an alternating electric field that separates the ions depending on their mobility, which is different for each type of molecule. Finally, to detect the presence of a specific disease biomarker, a 'compensation voltage' is applied that steers the molecules of interest to the detector, where they are counted.

The FAIMS chip is highly sensitive, able to detect at low parts per billion (ppb) and in some cases part per trillion (ppt) levels making it well suited for detecting breath VOCs which are typically present around the tens of ppb range. FAIMS is also highly selective, so even in a typical breath sample containing hundreds of VOCs unrelated to the disease of interest, it can ignore this chemical noise and detect with a low incidence of false positives. In addition the small size and software programmability of the FAIMS chip makes it suitable for near patient or point of care operation.



Early detection offers the best chance of survival in lung cancer



Lung cancer: enabling earlier intervention


These recent advances, both in breath sample collection and also in VOC analytical tools, open up new opportunities for widespread use of breath biopsy for early, non-invasive detection of disease, particularly for diseases that are currently poorly served. For example, lung cancer is the most common cancer in the world with an estimated 1.8 million cases and 1.6 million deaths in 2012¹. It is thought that in 2016 there will be more than 221,000 new cases of lung cancer in the United States, and with more than 158,000 deaths predicted, it is the number one cause of cancer deaths. Breast cancer, in comparison, has had a screening programme in place for some time, and the number of diagnoses is estimated to be 250,000, with but only 40,000 associated deaths².

The gold standard for lung cancer screening is low-dose computed tomography (LDCT), an imaging procedure that uses special low-dose x-ray equipment to create detailed scans of areas inside the body. In 2012, the National Lung Screening Trial (NLST) reported a mortality reduction of 20% compared to chest x-ray for lung cancer screening. The performance of the trial, however, has raised concerns regarding the high levels of false positives. For every 100,000 people screened for three years about 35,600 (35.6%) will have a false positive result, of these 1,800 would have an invasive procedure (5.1%) and 40 (1.1%) would have a major complication such as bleeding in the lung, a collapsed lung or an infection.

A breath biopsy test performed on patients who show positive for the LDCT test could greatly improve this; a breath test with 80% specificity would result in more than 13,500 fewer follow-up procedures for every 100,000 patients having a positive LDCT result. In addition, the fast turn around (days) of a breath test would greatly reduce the time patients would otherwise spend worrying about indeterminate LDCT test results, which under the current standards could take up to two years to obtain a definitive diagnosis. Alternatively, pre-screening the LDCT test population with a lung cancer breath biopsy test would similarly improve sensitivity and specificity while reducing patient discomfort.

The NHS and Innovate UK are currently funding the largest breath-based clinical trial in the world in order to investigate early stage lung cancer detection. Under chief investigator Robert Rintoul, the multi-centre Lung Cancer Indicator Detection (LuCID) study will sample up to 3,000 patients across the UK, Europe and the US referred under suspicion of lung cancer. Of this group 45% are expected to have lung cancer, about 7% of which are expected to be stage 1.

Besides early cancer detection, this study will also generate evidence to help guide downstream treatment decisions. Firstly, breath analysis could provide a non-invasive means to discern different tumour types, as the detected biomarkers strongly relate to the primary cell type and driving mutations. This could be particularly valuable in the



substantial fraction of patients with difficult to biopsy tumours. Secondly, breath biomarkers could potentially predict response to therapy, as has been suggested by more invasive metabolic profiling approaches. In that case ineffective treatment could be either prevented or discontinued quickly after identifying a lack of response. Finally, this approach could help to monitor patients after being treated with curative intent for a recurrence of their tumour.

This approach has benefits compared to CT-based follow-up as it can be done frequently without incurring significant healthcare costs or radiation exposure to patients. Taken together, non-invasive metabolomics through a breath biopsy has the potential to provide detailed insight into tumour biology impacting a wide range of aspects of diagnosis, treatment and monitoring.

Asthma: improving drug treatment regimes

The ability to measure disease activity by detecting downstream metabolites on the patient's breath has many applications to treatment stratification as well. For example, asthma is the most common chronic airway disease in the world. It is characterised by a combination of recurrent reversible air-flow limitation and characteristic symptoms as a consequence of chronic inflammation in the small airways. There are many different underlying mechanisms responsible for this inflammation, each requiring a different treatment. To date there is no test available that matches a patient to the treatment they respond to optimally other than a trial and error approach with progressive escalation towards higher dose treatment. Furthermore, up to one in five cases of asthma are classified as 'problematic' due to either symptom severity, control or both.

The difference in cell activity related to the different subtypes of asthma is also reflected by the metabolites these cells produce. As mentioned above, biomarkers such as fractional exhaled Nitric Oxide (FeNO) are being used to help guide treatment decisions: FeNO reflects a protective biochemical pathway of the lungs but lacks specificity as it is a single biomarker affected by many processes other than asthma. In spite of this it has been used to measure corticosteroid responsiveness in severe asthma patients, and has been shown to help assess treatment adherence for these patients, as FeNO levels are responsive to the correct application of steroid inhalers.

To address the 4-10% of asthma patients that are not able to achieve control over their symp-

toms due to poor responsiveness to current treatment options, pharmaceutical companies are currently undergoing clinical trials for novel biologic asthma therapies. To date finding the correct patient population for these new therapies has been extremely difficult, and currently available biomarkers such as FeNO, blood eosinophils and periostin have not adequately stratified these severe asthma patients.

Data from the U-BIOPRED study presented at the European Respiratory Society meeting in September 2015 showed that a FAIMS sensor was up to 63% more accurate in stratifying asthma patients in a head-to-head comparison with other gas sensors and classic laboratory techniques such as GC-MS, with an 81% accuracy in stratifying asthma patients responsive to anti-IgE treatment (XOLAIR®/Omalizumab). Being able to guide asthma patients to their most effective therapy will reduce patient discomfort and improve mortality, while reducing emergency care and hospitalisation costs.

The NHS and Innovate UK are currently funding the early stage clinical trial Stratification of Asthma Treatment by Breath Analysis (STRATA) to further explore the opportunity for breath biopsy among patients with problematic asthma. Being able to reliably match precision therapies to the correct patient population will drive further drug discovery. Embedding breath biopsy testing into early stage clinical trials would allow pharmaceutical companies to perform biomarker discovery early on in the drug development process. This gives them the opportunity to develop companion diagnostics that correctly target these therapies, thereby both improving treatment efficacy within the intended population while also creating an easily deployable test for identifying that population. Doing this maximises the chance of achieving regulatory approval for therapies that would otherwise fail due to health economic pressures.

Future possibilities for chemical biomarkers

Building on these advances in disease detection and stratification by therapeutic response, there are further possibilities for using breath biopsy for treatment monitoring and validation. Comparing breath tests with those performed pre-treatment could provide confirmation of the treatment's efficacy by detecting changes in disease-related biochemical activity, or even by detecting the metabolite signatures of the drug itself. Monitoring response in these activity signatures can help confirm the initial diagnosis. Also, tracking these

References

- 1** World Health Organization. International Agency for Research on Cancer. GLOBOCAN 2012: Estimated Cancer Incidence, Mortality and Prevalence Worldwide in 2012. Lung Cancer.
- 2** US SEER Cancer Statistics Factsheets: Lung and Bronchus Cancer. National Cancer Institute. Bethesda, MD.

changes over time could provide an early warning to changes in the underlying disease state and provide immediate feedback on treatment efficacy.

Miniaturised, yet highly sensitive and selective, breath analysers such as the FAIMS chip open up the possibility for future widespread adoption of breath biomarkers to not only aid in the development of precision therapies and define the relevant patient populations, but also to manage the effective application of those therapies. **DDW**

Billy Boyle is CEO of Owlstone Medical. An engineering graduate from Cambridge University, he is one of the original co-founders of Owlstone Inc, and in 2016 led the process to spin out Owlstone Medical Ltd where he is founding CEO.

Dr Mike Murphy is head of Business Development at Owlstone Medical. He received his PhD in semiconductor Physics from Cambridge University. Since then he has worked for a series of Cambridge-based start-ups focused on technology research and development, and is an alumni of the Techstars London accelerator programme.

ADVERTISEMENT INDEX

Advanced Cell Diagnostics	24,27,31	Clontech Laboratories, Inc	49	PerkinElmer, Inc	IFC
Agilent Technologies, Inc	4	EMD Millipore	43	SLAS	IBC
Analytik Jena AG	73	Eurofins Pharma Discovery Services	16	Source BioScience plc	6
BioTek Instruments, Inc	8	HTStec Ltd	21	Taconic Biosciences, Inc	62
BMG Labtech GmbH	52,61	IntelliCyt Corporation	3	ThermoFisher Scientific	32
Charles River Laboratories, Inc	50-51	Labcyte, Inc	70	Wagner Medizin und Pharmatechnik	57
CISBIO International SA	23,69	Pall ForteBio, LLC	15	Waters Corporation	OBC

Breath Biopsy® Products

Bring VOC Biomarker research to your own lab with our complete range of products for breath collection and analysis.

Breath collection



ReCIVA™ Breath Sampler

ReCIVA provides an out-of-the-box means to collect breath samples reliably and reproducibly.

- Collect different breath fractions (enriched bronchial or end-tidal)
- Pre-concentrate a large volume of breath onto sorbent tubes, enabling even VOCs present at low levels to be detected
- Developed in conjunction with a multi-disciplinary group of 100 breath analysis experts

owlstonemedical.com/reciva




CASPER Clean Air Supply

CASPER provides a consistent, portable supply of clean air, ideal for use with the ReCIVA™ Breath Sampler.

- Ensure breath samples are not contaminated by external VOCs
- Remove the need for a clean air supply line
- Portable - to wherever the patient is in clinic

owlstonemedical.com/casper

VOC analysis



Lonestar VOC Analyzer

The Lonestar VOC Analyzer is an easy to use instrument for the detection of VOC biomarkers in clinical samples such as breath and other biospecimens.

- Rapid, sensitive, and selective VOC analyzer
- Based on proven FAIMS technology
- Biomarker discovery and targeted analysis

owlstonemedical.com/lonestar

"We have found FAIMS technology able to non-invasively detect a broad range of diseases with high sensitivity and specificity."

Prof. J. Covington, Warwick University





Breath Biopsy[®] Biomarker Services

Discover novel non-invasive biomarkers for cancer, inflammatory and infectious disease with our comprehensive Breath Biomarker Services.

Overview



Owlstone Medical's Biomarker Services enable academics, clinicians and pharmaceutical companies to identify novel breath biomarkers for diagnostics as well as precision medicine applications including therapeutic response, patient stratification and outcome prediction.


We offer a complete set of services for analysis of volatile organic compound (VOC) biomarkers in breath and other biospecimens. Our services include everything from consultation on optimal study design to biomarker discovery and external blinded validation of the test in the context of the clinical use case.

Breath collection

We ensure that breath samples are collected reliably and reproducibly by installing the ReCIVA™ Breath Sampler at clinical sites, and providing training on optimal breath collection. Breath samples are pre-concentrated on stable sorbent tubes and shipped to Owlstone Medical's service laboratory for analysis.

Biomarker discovery

Owlstone Medical offers the largest laboratory for breath analysis. VOC Biomarkers are analyzed using our proven FAIMS technology. Our analysis experts develop robust classification algorithms from rich spectral FAIMS data. Our exceptional analytical platforms, rigorous quality standards and experienced staff deliver reliable, reproducible and rapid results.



In silico verification and validation

Classification algorithms are verified *in silico* using our proprietary database containing thousands of patient VOC profiles, to provide early evidence of likely clinical performance. Finally, an external blinded test validation is performed in patient groups representative of the intended use population.

Biomarker Services

Identify novel non-invasive
breath biomarkers

owlstonemedical.com/services

